

Datum: 2022-09-12 Dnr: RAÄ-2021-1228 ÖLM 2021.34 Handläggare: Kaj Thuresson

Instrument Report: µ- XRF, Microscopy, SEM/EDS

Sample Identification Code RAÄ Dnr RAÄ-2021-1228 Object no. OLM-39358-771 and OLM-39358-772 Sample Description of sample: Small dark chips/pieces with a powdery white substance on

surface.

Material: Leather

Point of analysis



Picture 1. The two samples with their own Eppendorf tube

Purpose

Investigation of the composition of the powdery substance on the leather.

Method

Sample preparation The whole sample was at first analysed with μ -XRF. A smaller portion of each sample was then taken with a scalpel and transferred to carbon tape, mounted on a SEM-stub, for investigation with optical light microscopy and SEM/EDS.

Instrument Parameters µXRF

	ſ	Χ	μXRF	Artax 800,	Mo X-ray	tube with	polycapillary	lens,	Bruker;	Berlin,	Germany
--	---	---	------	------------	----------	-----------	---------------	-------	---------	---------	---------

X single point analysis (spot size <100µm)

line scan (lateral resolution <100µm)

- elemental 2D mapping
- quantification, MQuant Calib, Bruker; Berlin, Germany
- quantification with standards

Riksantikvarieämbetet Artillerigatan 33 Box 1114 621 22 Visby Tel 08-5191 8000 E-post riksant@raa.se Hemsida www.raa.se Org.nr 202100-1090 Plusgiro 59994-4

Bankgiro 5052-3620



Voltage 50 kV
Current 600 μA
Scan time per point 10 s
Number of measurement points 8 + 8
Filter 🗙 no filter 🗌 Al 315 μm 🗌 Mo 12.50 μm 🗌 other
Lens 0.060
Atmosphere X air He for light element detection

Qualitative Results µ-XRF



Picture 2. Sample OLM-39358-771 was analysed at eight µ-points. All spectra were comparable (no difference between light or dark areas on sample). Mostly calcium (Ca), small amounts of iron (Fe) and manganese (Mn), and possibly phosphor (P), were detected in addition to trace amounts of copper (Cu), zinc (Zn) and sulfur (S).



Picture 3. Sample OLM-39358-772 was analysed at eight μ-points. All spectra were comparable (no difference between light or dark areas on sample). Mostly iron (Fe) and calcium (Ca), small amounts of manganese (Mn), copper (Cu), potassium (K) and titanium (Ti) were detected, in addition to trace amounts of zinc (Zn), sulfur (S) and phosphor (P).



Microscopy results Microscope used: LEICA S9i External light source: VWR VisiLight[™] CL-150



Picture 4. OLM-39358-771 sample on carbon tape in two different magnifications.



Picture 5. OLM-39358-772 sample on carbon tape in two different magnifications



SEM/EDS results Instrument used: JEOL JSM-IT500



____ 20mm

Sample OLM-39358-771 and OLM-39358-772 in SEM chamber



Items	Value
measurement conditions	
Acceleration voltage	15.00 kV
Probe current	-
Magnification	x 600
Process time	T2
Measurement detector	First
Live time	58.85 seconds
Real time	60.09 seconds
Dead time	2.00
Count rate	4713.00 CPS

Signal BED-S Landing Voltage 15.0 kV WD 10.9 mm Magnification x600

Vacuum Mode Low Vacuum





■ 20µm





Items	Value	
measurement conditions		
Acceleration voltage	15.00 kV	
Probe current	-	
Magnification	x 2000	
Process time	T2	
Measurement detector	First	
Live time	38.20 seconds	
Real time	39.03 seconds	
Dead time	2.00	
Count rate	4846.00 CPS	

= 10 μm

Signal BED-S Landing Voltage 15.0 kV WD 10.9 mm Magnification x600 Vacuum Mode Low Vacuum

Element	Line	Mass%	Atom%
0	К	47.87±0.29	69.54±0.42
Р	К	1.40±0.05	1.05±0.03
Ca	К	50.73±0.27	29.41±0.15
Total		100.00	100.00
Spc_002			Fitting ratio 0.5913



SWEDISH NATIONAL HERITAGE BOARD RIKSANTIKVARIEÄMBETET



Signal BED-S Landing Voltage 15.0 kV WD 14.0 mm Magnification x30 F.O.V. 4.267 x 3.200 mm Probe Current Std. 65.0 Scan Rotation 191.4° Pressure 36 Pa



Signal BED-S

WD 13.8 mm

Magnification x600

Landing Voltage 15.0 kV

Vacuum Mode Low Vacuum

Vacuum Mode Low Vacuum

500 μm

Items	Value
measurement conditions	
Acceleration voltage	15.00 kV
Probe current	-
Magnification	x 600
Process time	T2
Measurement detector	First
Live time	43.92 seconds
Real time	44.78 seconds
Dead time	2.00
Count rate	4610.00 CPS

120 µm

Element	Line	Mass%	Atom%
С	К	52.81±0.10	60.69±0.12
0	К	43.95±0.19	37.92±0.17
Si	К	1.89±0.04	0.93±0.02
Ca	К	1.36±0.04	0.47±0.02
Total		100.00	100.00
Spc_005			Fitting ratio 0.5666

